



**Abbott**

# Quantisal™ Oral Fluid Testing

## IMPROVING YOUR SCREENING OPTIONS

We are pleased to announce we are advancing our suite of drug testing solutions by transitioning from OraSure Intercept® to the Abbott Quantisal™ Oral Fluid Testing device.

The Quantisal™ Oral Fluid Testing device provides the following benefits:

- The oral fluid immunoassays and cutoffs used are aligned to SAMHSA.
- Quantisal contains no artificial stimulants which could interfere with the test results.
- A unique volume adequacy indicator allowing for a measured sample collection.
- The Quantisal device supports your current drug menu and provides additional test options as shown below.

## KEY DATES

- The final date to order OraSure Intercept devices is 3/31/2024.
- The final date that OraSure Intercept specimens will be accepted for testing is 5/31/2024.
- The chain of custody forms that you use today for OraSure Intercept will still be accepted with an Abbott Quantisal specimen. Any chain of custody forms ordered after March 31, 2024 will reflect your new Abbott Quantisal panel.

**The oral fluid immunoassays used with Quantisal are developed and launched with cutoffs aligned to SAMHSA.**

LAB DRUG MENU	CONFIRMATION MENU	SCREEN LEVEL NG/ML	RECOMMENDED CONFIRM LEVEL NG/ML
Amphetamine	Amphetamine, MDA	50	25
Benzodiazepines	Benzodiazepines (Alprazolam, Clonazepam, Diazepam, Lorazepam, Nordiazepam, Oxazepam, Temazepam)	20	5
Barbiturates	Barbiturates (Amobarbital, Butalbital, Pentobarbital, Phenobarbital, Secobarbital)	50	50
Buprenorphine	Buprenorphine, Norbuprenorphine	5	3
Cocaine	Cocaine, Benzoylcegonine	15	8
Cotinine	Cotinine	50	30
Ethanol (Alcohol)	Ethanol	0.04	0.02
Marijuana	THC (Δ9-Tetrahydrocannabinol)	4	2
Methadone	Methadone, EDDP	50	15
Methamphetamine	Methamphetamine, MDMA	50	25
Opiates	Codeine, Morphine	30	15
Extended Opiates	Codeine, Morphine, Hydrocodone, Hydromorphone	30	15
Phencyclidine, PCP	Phencyclidine, PCP	10	10
<b>NOW AVAILABLE AS INDIVIDUAL ANALYTES</b>			
Oxycodone/Oxymorphone*	Oxycodone, Oxymorphone	30	15
6AM - 6-Acetylmorphine*	6-Acetylmorphine	4	2

\*Oxycodone and 6AM are now available as individual tests to enhance your drug testing menu with Oral Fluid

**PLEASE CONTACT CUSTOMER SERVICE IF YOU HAVE ANY QUESTIONS OR WANT ADDITIONAL INFORMATION REGARDING THIS NEW PRODUCT TRANSITION.**